

Product Datasheet

GFAP (labeling) antibody (orb1471989)

Description	Mouse Mouse antibody to GFAP (labeling).
Species/Host	Mouse
Conjugation	Unconjugated
Target	GFAP
Preservatives	PBS pH 7.4
Form/Appearance	Liquid
Storage	-20°C
Note	For research use only
Application notes	Suggested applications: EIA, CLIA, RIA, IHC, WB
Clonality	Monoclonal
Source	Ascites
Expiration Date	12 months from date of receipt.