

Product Datasheet

BACE1 (labeling) antibody (orb1471993)

Description	Mouse Mouse antibody to BACE1 (labeling).
Species/Host	Mouse
Conjugation	Unconjugated
Target	BACE1
Preservatives	PBS pH 7.4
Form/Appearance	Liquid
Storage	-20°C
Note	For research use only
Application notes	Suggested applications: EIA, CLIA, RIA, IHC, WB
Clonality	Monoclonal
Source	Ascites
Expiration Date	12 months from date of receipt.