

Product Datasheet

FKBP52 antibody (orb147210)

Description

Mouse monoclonal to FKBP52 (APC). Hsp90 is

the fo...

Species/Host

Mouse

Reactivity

Canine, Hamster, Human, Mouse, Rat

Conjugation

APC

Tested

ICC, IHC, IP, WB

Applications

Immunogen

Synthetic peptide corresponding to the residues of human FKBP52

Target

FKBP52

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

0.5 µg/ml was sufficient for detection of FKBP52 in 20 µg total protein using WB by colorimetric immunoblot analysis using Goat Anti-Mouse IgG:HRP as the secondary.

Isotype

IgG2a

Clonality

Monoclonal

MW

52kDa

Uniprot ID

[Q02790](#)

NCBI

[NP_002005.1](#)

Entrez

[2288](#)

Dilution Range

WB (1:2000), IHC (1:250), ICC/IF (1:1000), IP (5µg)

Expiration Date

12 months from date of receipt.