



Product Datasheet

LAMP1 antibody [Out of stock] (orb147219)

Description	Mouse monoclonal to LAMP1 (ATTO-488). Lysosome associated membrane proteins, or LAMP1 and LAMP2, are ...
Species/Host	Mouse
Reactivity	Hamster, Human, Mouse, Rat
Conjugation	ATTO-488
Tested Applications	ICC, IP, WB
Immunogen	Rat liver lysosomal membrane preparations
Target	LAMP1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml was sufficient for detection of LAMP1 in rat liver microsome by ECL immunoblot analysis.
Isotype	IgG1
Clonality	Monoclonal
MW	120kDa
Uniprot ID	P14562
NCBI	NP_036989.1
Entrez	25328
Dilution Range	WB (1:1000), ICC/IF (1:1000)
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of human
tra...



Western Blot
analysis of Rat liver
micro...