

Product Datasheet

LAMP1 antibody [Out of stock] (orb147221)

Description

Mouse monoclonal to LAMP1 (ATTO-594).

or LAMP1 and LAMP2, are ...

Species/Host

Mouse

Reactivity

Hamster, Human, Mouse, Rat

Conjugation

ATTO-594

Tested

ICC, IP, WB

Applications
Immunogen

Rat liver lysosomal membrane preparations

Target

LAMP1

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml was sufficient for detection of LAMP1 in rat liver microsome by ECL immunoblot analysis.

Isotype

IgG1

Clonality

Monoclonal

MW

120kDa

Uniprot ID

[P14562](#)

NCBI

[NP_036989.1](#)

Entrez

[25328](#)

Dilution Range

WB (1:1000), ICC/IF (1:1000)

Expiration Date

12 months from date of receipt.