

Product Datasheet

LAMP1 antibody (orb147227)

Description	Mouse monoclonal to LAMP1 (APC). Lysosome associated membrane proteins, or LAMP1 and LAMP2, are major...
Species/Host	Mouse
Reactivity	Hamster, Human, Mouse, Rat
Conjugation	APC
Tested Applications	ICC, IP, WB
Immunogen	Rat liver lysosomal membrane preparations
Target	LAMP1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml was sufficient for detection of LAMP1 in rat liver microsome by ECL immunoblot analysis.
Isotype	IgG1
Clonality	Monoclonal
MW	120kDa
Uniprot ID	P14562
NCBI	NP_036989.1
Entrez	25328
Dilution Range	WB (1:1000), ICC/IF (1:1000)
Expiration Date	12 months from date of receipt.