

Product Datasheet

Alpha A Crystallin antibody [Out of stock] (orb147253)

Description

Mouse monoclonal to Alpha A

	crystallins are major water-soluble len...
Species/Host	Mouse
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	ATTO-488
Tested Applications	ELISA, ICC, WB
Immunogen	Native Alpha Crystallin
Target	Alpha A Crystallin
Preservatives	PBS pH 7.2, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	0.5 µg/ml was sufficient for detection of 100 ng purified alphaA crystalline by colorimetric immunoblot analysis using Goat Anti-Mouse IgG:HRP as the secondary.
Isotype	IgG1
Clonality	Monoclonal
MW	20kDa
Uniprot ID	P02489
NCBI	NP_000385.1
Entrez	1409
Dilution Range	WB (1:2000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.