

Product Datasheet

Phospho-RSK1 p90 (Thr359/Ser363) Rabbit mAb Antibody (orb1472614)

Description	Phospho-RSK1 p90 (Thr359/Ser363) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q15418
Expiration Date	12 months from date of receipt.