

Product Datasheet

HSP40, Hdj1 antibody [Out of stock] (orb147270)

Description

Mouse monoclonal to Hsp40, Hdj1

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-488

Tested

ELISA, ICC, IHC, IP, WB

Applications
Immunogen

Recombinant Protein HSP40 (Hdj1)

Target

HSP40

Preservatives

PBS pH 7.2, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

0.5 µg/ml of SMC-145 was sufficient for detection of HSP40 (HDJ1) in 15 µg of HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

40kDa

Uniprot ID

[P25685](#)

NCBI

[NP_006136.1](#)

Entrez

[3337](#)

Dilution Range

WB (1:2000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.