

Product Datasheet

HSP40, YDJ1 antibody [Out of stock] (orb147321)

Description

Mouse monoclonal to Hsp40, YDJ1

proteins comprise a large protein
famil...

Species/Host

Mouse

Reactivity

Yeast

Conjugation

ATTO-488

Tested

ELISA, IP, WB

Applications
Immunogen

Full length protein yeast HSP40 (YDJ1)

Target

HSP40

Preservatives

50% glycerol, 0.09% sodium azide
*Storage buffer may change when
conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored
at 4°C

Note

For research use only

Application notes

0.5 µg/ml of SMC-150 was sufficient for
detection of 50 ng YDJ1 by colorimetric
immunoblot analysis using Goat anti-
mouse IgG:HRP as the secondary
antibody.

Isotype

IgG1 Kappa

Clonality

Monoclonal

MW

40kDa

Uniprot ID
[P25491](#)
NCBI
[NP_014335.1](#)
Entrez
[855661](#)
Dilution Range

WB (1:2000)

Expiration Date

12 months from date of receipt.