

Product Datasheet

Myc-tag Blocking Peptide (orb1473254)

Description	Myc-tag Blocking Peptide
Tested Applications	Blocking
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Expiration Date	6 months from date of receipt.