


Product Datasheet

Swine IL-6 ELISA Kit (orb1473286)

Description	Swine IL-6 ELISA Kit
Reactivity	Porcine
Range	125 pg/ml - 8000 pg/ml
Tested Applications	ELISA
Target	IL6
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P26893
Sensitivity	60 pg/ml
Entrez	399500
Expiration Date	Please enquire.



Representative
Standard
Curve