


Product Datasheet

Swine IL-12/IL-23p40 ELISA Kit (orb1473289)

Description	Swine IL-12/IL-23p40 ELISA Kit
Reactivity	Porcine
Range	78 pg/ml - 5000 pg/ml
Tested Applications	ELISA
Target	IL12A; IL12B
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q28938, Q29053
Sensitivity	40 pg/ml
Entrez	397076, 397053
Expiration Date	Please enquire.



Representative
Standard
Curve