


Product Datasheet

Swine CRP ELISA Kit (orb1473290)

Description	Swine CRP ELISA Kit
Reactivity	Porcine
Range	78 pg/ml - 5000 pg/ml
Tested Applications	ELISA
Target	CRP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	O19062
Sensitivity	40 pg/ml
Entrez	396842
Expiration Date	Please enquire.



Representative
Standard
Curve