


Product Datasheet

Rat TRK B ELISA Kit (orb1473294)

Description	Rat TRK B ELISA Kit
Reactivity	Rat
Range	93.6 pg/ml - 6000 pg/ml
Tested Applications	ELISA
Target	NTRK2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Cell Lysates
Sample Types	Cell Culture Supernates, Cell Lysates
Uniprot ID	Q63604
Sensitivity	50 pg/ml
Entrez	25054
Expiration Date	Please enquire.



Representative
Standard
Curve