

Product Datasheet

HSP40, YDJ1 antibody (orb147330)

Description

Mouse monoclonal to Hsp40, YDJ1 (Biotin).

protein family,...

Species/Host

Mouse

Reactivity

Yeast

Conjugation

Biotin

Tested

ELISA, IP, WB

Applications
Immunogen

Full length protein yeast HSP40 (YDJ1)

Target

HSP40

Preservatives

50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

0.5 µg/ml of SMC-150 was sufficient for detection of 50 ng YDJ1 by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1 Kappa

Clonality

Monoclonal

MW

40kDa

Uniprot ID

[P25491](#)

NCBI

[NP_014335.1](#)

Entrez

[855661](#)

Dilution Range

WB (1:2000)

Expiration Date

12 months from date of receipt.