


Product Datasheet

Rat Plexin A1 ELISA Kit (orb1473303)

Description	Rat Plexin A1 ELISA Kit
Reactivity	Rat
Range	312 pg/ml - 20000 pg/ml
Tested Applications	ELISA
Target	PLXNA1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	D3Z981
Sensitivity	150 pg/ml
Entrez	362398
Expiration Date	Please enquire.



Representative
Standard
Curve