


Product Datasheet

Rat FGFBP1 ELISA Kit (orb1473316)

Description	Rat FGFBP1 ELISA Kit
Reactivity	Rat
Range	78 pg/ml - 5000 pg/ml
Tested Applications	ELISA
Target	FGFBP1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9QY10
Sensitivity	40 pg/ml
Entrez	64535
Expiration Date	Please enquire.



Representative
Standard
Curve