


Product Datasheet

Mouse TIM3 ELISA Kit (orb1473338)

Description	Mouse TIM3 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	HAVCR2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q8VIM0
Sensitivity	15 pg/ml
Entrez	171285
Expiration Date	Please enquire.



Representative
Standard
Curve