


Product Datasheet

Mouse SMOC1 ELISA Kit (orb1473340)

Description	Mouse SMOC1 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	SMOC1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q8BLY1
Sensitivity	15 pg/ml
Entrez	64075
Expiration Date	Please enquire.



Representative
Standard
Curve