


Product Datasheet

Mouse Serpin E2 ELISA Kit (orb1473343)

Description	Mouse Serpin E2 ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20000 pg/ml
Tested Applications	ELISA
Target	SERPINE2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q07235
Sensitivity	150 pg/ml
Entrez	20720
Expiration Date	Please enquire.



Representative
Standard
Curve