


Product Datasheet

Mouse Serpin D1 ELISA Kit (orb1473344)

Description	Mouse Serpin D1 ELISA Kit
Reactivity	Mouse
Range	780 pg/ml - 50000 pg/ml
Tested Applications	ELISA
Target	SERPIND1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P49182
Sensitivity	400 pg/ml
Entrez	15160
Expiration Date	Please enquire.



Representative
Standard
Curve