


Product Datasheet

Mouse Serpin A1b ELISA Kit (orb1473345)

Description	Mouse Serpin A1b ELISA Kit
Reactivity	Mouse
Range	1560 pg/ml - 100000 pg/ml
Tested Applications	ELISA
Target	SERPINA1B
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P22599
Sensitivity	800 pg/ml
Entrez	20701
Expiration Date	Please enquire.



Representative
Standard
Curve