


Product Datasheet

Mouse LIF ELISA Kit (orb1473359)

Description	Mouse LIF ELISA Kit
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	LIF
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P09056
Sensitivity	4 pg/ml
Entrez	16878
Expiration Date	Please enquire.



Representative
Standard
Curve