


Product Datasheet

Mouse HRG ELISA Kit (orb1473365)

Description	Mouse HRG ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	HRG
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9ESB3
Sensitivity	80 pg/ml
Entrez	94175
Expiration Date	Please enquire.



Representative
Standard
Curve