


Product Datasheet

Mouse GDF8 ELISA Kit (orb1473369)

Description	Mouse GDF8 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	MSTN
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	008689
Sensitivity	15 pg/ml
Entrez	17700
Expiration Date	Please enquire.



Representative
Standard
Curve