


Product Datasheet

Mouse Gc-globulin ELISA Kit (orb1473370)

Description	Mouse Gc-globulin ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20000 pg/ml
Tested Applications	ELISA
Target	GC
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma, Cell Lysates
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates
Uniprot ID	P21614
Sensitivity	150 pg/ml
Entrez	14473
Expiration Date	Please enquire.



Representative
Standard
Curve