


Product Datasheet

Mouse CXCL5 ELISA Kit (orb1473377)

Description	Mouse CXCL5 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	CXCL5
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates, Urine
Uniprot ID	P50228
Sensitivity	8 pg/ml
Entrez	20311
Expiration Date	Please enquire.



Representative
Standard
Curve