


Product Datasheet

Mouse Complement C1R ELISA Kit (orb1473381)

Description	Mouse Complement C1R ELISA Kit
Reactivity	Mouse
Range	625 pg/ml - 40000 pg/ml
Tested Applications	ELISA
Target	C1RA; C1RB
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q8CG16 , Q8CFG9
Sensitivity	300 pg/ml
Entrez	667277 , 50909
Expiration Date	Please enquire.



Representative
Standard
Curve