


Product Datasheet

Mouse CD307e ELISA Kit (orb1473384)

Description	Mouse CD307e ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	FCRL5
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q68SN8
Sensitivity	15 pg/ml
Entrez	329693
Expiration Date	Please enquire.



Representative
Standard
Curve