


Product Datasheet

Mouse BR Cadherin ELISA Kit (orb1473391)

Description	Mouse BR Cadherin ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	CDH12
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q5RJH3
Sensitivity	80 pg/ml
Entrez	215654
Expiration Date	Please enquire.



Representative
Standard
Curve