


Product Datasheet

Mouse BMP9 ELISA Kit (orb1473392)

Description	Mouse BMP9 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	GDF2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Bone tissue, Cell Culture Supernates, Serum, Plasma
Sample Types	Bone tissue, Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9WV56
Sensitivity	8 pg/ml
Entrez	12165
Expiration Date	Please enquire.



Representative
Standard
Curve