


Product Datasheet

Mouse BMP6 ELISA Kit (orb1473393)

Description	Mouse BMP6 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	BMP6
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Bone tissue, Cell Culture Supernates, Serum
Sample Types	Bone tissue, Cell Culture Supernates, Serum
Uniprot ID	P20722
Sensitivity	80 pg/ml
Entrez	12161
Expiration Date	Please enquire.



Representative
Standard
Curve