


Product Datasheet

Human Neudesin ELISA Kit (orb1473441)

Description	Human Neudesin ELISA Kit
Reactivity	Human
Range	93.6 pg/ml - 6000 pg/ml
Tested Applications	ELISA
Target	NENF
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9UMX5
Sensitivity	50 pg/ml
Entrez	29937
Expiration Date	Please enquire.



Representative
Standard
Curve