


Product Datasheet

Human HSPG2 ELISA Kit (orb1473459)

Description	Human HSPG2 ELISA Kit
Reactivity	Human
Range	625 pg/ml - 40000 pg/ml
Tested Applications	ELISA
Target	HSPG2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P98160
Sensitivity	300 pg/ml
Entrez	3339
Expiration Date	Please enquire.



Representative
Standard
Curve