


Product Datasheet

Human HO-1 ELISA Kit (orb1473461)

Description	Human HO-1 ELISA Kit
Reactivity	Human
Range	93.6 pg/ml - 6000 pg/ml
Tested Applications	ELISA
Target	HMOX1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	P09601
Sensitivity	50 pg/ml
Entrez	3162
Expiration Date	Please enquire.



Representative
Standard
Curve