


Product Datasheet

Human GDF8 ELISA Kit (orb1473466)

Description	Human GDF8 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	MSTN
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	O14793
Sensitivity	8 pg/ml
Entrez	2660
Expiration Date	Please enquire.



Representative
Standard
Curve