


Product Datasheet

Human Complement C1S ELISA Kit (orb1473491)

Description	Human Complement C1S ELISA Kit
Reactivity	Human
Range	780 pg/ml - 50000 pg/ml
Tested Applications	ELISA
Target	C1S
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P09871
Sensitivity	400 pg/ml
Entrez	716
Expiration Date	Please enquire.



Representative
Standard
Curve