


Product Datasheet

Human CNDP1 ELISA Kit (orb1473495)

Description	Human CNDP1 ELISA Kit
Reactivity	Human
Range	3125 pg/ml - 200000 pg/ml
Tested Applications	ELISA
Target	CNDP1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q96KN2
Sensitivity	1500 pg/ml
Entrez	84735
Expiration Date	Please enquire.



Representative
Standard
Curve