


Product Datasheet

Human CD120a ELISA Kit (orb1473523)

Description	Human CD120a ELISA Kit
Reactivity	Human
Range	12.5 pg/ml - 800 pg/ml
Tested Applications	ELISA
Target	TNFRSF1A
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P19438
Sensitivity	7 pg/ml
Entrez	7132
Expiration Date	Please enquire.



Representative
Standard
Curve