


Product Datasheet

Human Cathepsin L ELISA Kit (orb1473529)

Description	Human Cathepsin L ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	CTSL
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma, Cell Lysates
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates
Uniprot ID	P07711
Sensitivity	30 pg/ml
Entrez	1514
Expiration Date	Please enquire.



Representative
Standard
Curve