

Product Datasheet

Mouse Human IgA antibody (FITC) (orb1473579)

Description	Mouse polyclonal antibody against Mouse Human IgA conjugated to
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	ELISA, FC, IF, IH
Immunogen	Human IgA
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Mouse
Dilution Range	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000)
Expiration Date	12 months from date of receipt.