

Product Datasheet

Goat Human IgG+IgM+IgA antibody (Biotin) (orb1473594)

Description	Goat polyclonal antibody against Goat Human IgG+IgM+IgA conjugated to
Species/Host	Goat
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, IH, WB
Immunogen	Human IgG+IgM+IgA
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Goat
Dilution Range	E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.