


Product Datasheet

TCF1 alpha antibody (orb1473638)

Description	Mouse monoclonal antibody to TCF1 alpha
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TCF1 alpha. The exact sequence is proprietary.
Target	LEF1
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q9UJU2
Entrez	51176
Dilution Range	WB (1/1000 - 1/2000)
Expiration Date	12 months from date of receipt.



Western blot analysis of TCF1 alpha expr...