

Product Datasheet

SPAK antibody (orb1473643)

Description	Mouse monoclonal antibody to SPAK
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SPAK. The exact sequence is proprietary.
Target	STK39
Preservatives	Mouse IgG1. Supplied in crude ascites with 0.01% sodium azide.
Form/Appearance	Mouse IgG1. Supplied in crude ascites with 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q9UEW8
Entrez	27347
Dilution Range	WB (1/500 - 1/2000)
Expiration Date	12 months from date of receipt.



Western blot analysis of SPAK expression...