

Product Datasheet

FN3K antibody (orb1473676)

Description	Mouse monoclonal antibody to FN3K
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human FN3K. The exact sequence is proprietary.
Target	FN3K
Preservatives	Mouse IgG1. Supplied in crude ascites with 0.01% sodium azide.
Form/Appearance	Mouse IgG1. Supplied in crude ascites with 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q9H479
Entrez	64122
Dilution Range	WB (1/500 - 1/4000)
Expiration Date	12 months from date of receipt.

