

Product Datasheet

UBL5 antibody (orb1473684)

Description	Mouse monoclonal antibody to UBL5
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human UBL5. The exact sequence is proprietary.
Target	UBL5
Preservatives	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q9BZL1 , Q9EPV8
Entrez	59286 , 66177
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.

