

Product Datasheet

p38 alpha (MAP Kinase) antibody (orb147369)

Description

Mouse monoclonal to p38 alpha (RPE).

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

RPE

Tested

AM, ELISA, IHC, IP, WB

Applications
Immunogen

Full length recombinant protein expressed in E.coli cells

Target

p38 MAPK

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

Detects ~38kDa protein corresponding to p38α MAPK when loaded with 6 ng of purified p38α by chemiluminescent immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

38kDa

Uniprot ID

[Q16539](#)

NCBI

[NP_001306.1](#)

Entrez

[1432](#)

Dilution Range

WB (1:1000), IHC (1:1000)

Expiration Date

12 months from date of receipt.