



Product Datasheet

CRMP-5 antibody (orb1473693)

Description	Mouse monoclonal antibody to CRMP-5
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human CRMP-5. The exact sequence is proprietary.
Target	DPYSL5
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q9JHU0 , Q9EQF6 , Q9BPU6
Entrez	65208 , 100910964 , 65254 , 56896
Dilution Range	WB (1/1000 - 1/2000), IF/IC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of CRMP-5 expression...



Immunofluorescent analysis of CRMP-5 staining...