


Product Datasheet

PURB antibody (orb1473704)

Description	Mouse monoclonal antibody to PURB
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PURB. The exact sequence is proprietary.
Target	PURB
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	O35295 , Q96QR8
Entrez	19291 , 5814
Dilution Range	WB (1/500 - 1/4000)
Expiration Date	12 months from date of receipt.



Western blot analysis of PURB expression...